

# PIILOSET PETSET

*Preservative-free multidose comfort eye drop*

Dry eyes syndrome (Keratoconjunctivitis sicca) is rather a common disorder among humans and animals. It can be caused by, for instance, lacrimal gland dysfunction, infections, certain medications, neurologic diseases etc.

PetSet Comfort drops relieve unpleasant dry eye symptoms easily and naturally. PetSet protects and moisturizes the eye dryness, itching the eye and decreases irritation and discomfort.

## NATURAL RELIEF

The active ingredient of PetSet Comfort drop is 0.2% hyaluronan. Hyaluronan is an eye-friendly moisturizer with muco-mimetic and muco-adhesive properties that improve tear film stability and keep the ocular surface comfortable. To ensure a natural feeling in the eyes, chemical and physical properties of PetSet Comfort drop are adjusted to resemble the healthy tear film.

## GENTLE

Preservative-free PetSet Comfort drop suits also sensitive eyes and prolonged use. Free from preservatives and ingredients of animal or plant origin. PetSet Comfort drop has not undergone animal testing.

## PACKING

PetSet comfort drop is packed into a 10 ml (approx 200 drops) multidose dispensing system designed for sterile and preservative-free ophthalmic products. The dispenser is calibrated for the properties of PetSet comfort drop. The dispensing system keeps the product sterile for 3 months from opening.

## APPLICATION

Apply 1-2 drops into the eye when needed.



10 ml multidose dispenser  
Shelf-life 2 years



Finnsusp Ltd | Pääskykalliontie 5, 21420 Lieto, Finland | Tel +358 2 489 110 | info@piiloseet.fi | www.piiloseet.fi



## PETSET EYE BATH

*For safe cleansing of eyes and eyelids*

PetSet is a gentle eye bath for animals. The sterile solution includes eye friendly cleansing ingredients that effectively remove impurities from the eyes and eyelids. Fragrance-free PetSet does not irritate eyes. Furthermore, preservatives are set to optimally low level.

### INGREDIENTS

Sodium Borate, Sodium chloride, Hydrochloric acid, DiSodium EDTA, EDG, Biguanide.

PetSet has not been tested on animals and does not contain ingredients of animal origin. Developed and manufactured in Finland.

### INSTRUCTIONS

Apply a few drops of PetSet gently into the eye and let the dog/cat blink.

Gently wipe the excess bath with a clean cotton pad around the eye.

Suitable for daily use.

### AVAILABILITY

Bottle size 120 ml

Discard after 4 months from opening

EXP 3 years



Gentle Care  
for Your Pet



No irritation!

- Fragrance-free
- Preservative set on optimally low level



Finnsusp Ltd | Pääskykalliontie 5, 21420 Lieto, Finland | Tel +358 2 489 110 | info@piilose.fi | www.piilose.fi

